

Product datasheet for **TP725206**

TEV Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant TEV Protease
Species:	TEV
Expression cDNA Clone or AA Sequence:	Gly2038-Asn2273
Tag:	N-6His
Buffer:	Supplied as a 0.2 um filtered solution of 25mMTris, 0.5mM EDTA, 2.5mM DTT, 25% glycerol, 0.05% TritonX-100, pH7.5
Note:	Recombinant TEV Nuclear Inclusion protein A is produced by our E.coli expression system and the target gene encoding Gly2038-Asn2273 is expressed with a 6His tag at the N-terminus.
Stability:	12 months from date of despatch
Summary:	TEV Nuclear Inclusion protein A has RNA-binding and proteolytic activities. It can form a dimer with disulfide-link.



[View online »](#)