

Product datasheet for **TP725208**

RAB23 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Ras-related protein Rab-23
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-G203
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Ras-related protein Rab-23 is produced by E.coli expression system and the target gene encoding M1-G203 is expressed.
Stability:	12 months from date of despatch
Locus ID:	51715
UniProt ID:	Q9ULC3



[View online »](#)