

Product datasheet for **TP725209**

UPF3A Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human UPF3A Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	E327-E476
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human UPF3A Protein is produced by E.coli expression system and the target gene encoding E327-E476 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	65110
UniProt ID:	Q9H1J1
Synonyms:	Regulator of nonsense transcripts 3A; Nonsense mRNA reducing factor 3A; Up-frameshift suppressor 3 homolog A; hUpf3; RENT3A; UPF3



[View online »](#)