

Product datasheet for **TP725213**

MRRF Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human RRF
Species:	Human
Expression cDNA Clone or AA Sequence:	A56-G262
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human RRF is produced by E.coli expression system and the target gene encoding A56-G262 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	92399
UniProt ID:	Q96E11
Synonyms:	Ribosome-recycling factor, mitochondrial; Ribosome-releasing factor, mitochondrial; MRRF
Protein Families:	Transmembrane



[View online »](#)