

Product datasheet for **TP725214**

RBFOX1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human RBFOX1 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-F128
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human RBFOX1 Protein is produced by E.coli expression system and the target gene encoding M1-F128 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	54715
UniProt ID:	<u>Q9NWB1</u>
Synonyms:	RNA binding protein fox-1 homolog 1; Ataxin-2-binding protein 1; Fox-1 homolog A; Hexaribonucleotide-binding protein 1; A2BP; A2BP1; FOX1; HRNBP1



[View online »](#)