

Product datasheet for **TP725224**

SESN1 Human Recombinant Protein

Product data:

| | |
|---------------------------------------|--|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant Human Sestrin-1/SESN1(N-His-SUMO) |
| Species: | Human |
| Expression cDNA Clone or AA Sequence: | T101-V200 |
| Tag: | N-His-SUMO |
| Buffer: | Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4 |
| Note: | MADE TO ORDER Recombinant Human Sestrin-1/SESN1 is produced by E.coli expression system and the target gene encoding T101-V200 is expressed with a His&SUMO tag at the N-terminus. |
| Stability: | 12 months from date of despatch |
| Locus ID: | 27244 |
| UniProt ID: | Q9Y6P5 |



[View online »](#)