

Product datasheet for **TP725225**

STAC Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human STAC Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	L209-S318
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human STAC Protein is produced by E.coli expression system and the target gene encoding L209-S318 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	6769
UniProt ID:	Q99469
Synonyms:	SH3 and cysteine-rich domain-containing protein; Src homology 3 and cysteine-rich domain-containing protein; STAC1



[View online »](#)