

Product datasheet for **TP725226**

SIPA1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Sipa-1(N-His-SUMO)
Species:	Human
Expression cDNA Clone or AA Sequence:	E949-A1042
Tag:	N-His-SUMO
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Sipa-1 is produced by E.coli expression system and the target gene encoding E949-A1042 is expressed with a His&SUMO tag at the N-terminus.
Stability:	12 months from date of despatch
Locus ID:	6494
UniProt ID:	Q96FS4



[View online »](#)