

Product datasheet for **TP725233**

DENTT (TSPYL2) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human TSPY-like protein 2
Species:	Human
Expression cDNA Clone or AA Sequence:	W344-G693
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human TSPY-like protein 2 is produced by E.coli expression system and the target gene encoding W344-G693 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	64061
UniProt ID:	Q9H2G4
Synonyms:	Testis-specific Y-encoded-like protein 2; Cell division autoantigen 1; Cutaneous T-cell lymphoma-associated antigen se20-4;Differentially-expressed nucleolar TGF-beta1 target protein; CTCL-associated antigen se20-4; Nuclear protein of 79 kDa ; NP79; TSPYL2; CDA1; DENTT; TSPX
Protein Families:	Druggable Genome



[View online »](#)