

## Product datasheet for **TP725234**

### THRSP Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human THRSP (N-His-SUMO)
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-W146
Tag:	N-His-SUMO
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human THRSP is produced by E.coli expression system and the target gene encoding M1-W146 is expressed with a His&SUMO tag at the N-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	7069
UniProt ID:	<a href="#">Q92748</a>
Synonyms:	Thyroid hormone-inducible hepatic protein; Spot 14 protein; S14; SPOT14
Protein Families:	Transcription Factors



[View online »](#)