

Product datasheet for **TP725236**

TRIT1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human tRNA dimethylallyltransferase/TRIT1
Species:	Human
Expression cDNA Clone or AA Sequence:	E357-H446
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human tRNA dimethylallyltransferase/TRIT1 is produced by E.coli expression system and the target gene encoding E357-H446 is expressed.
Stability:	12 months from date of despatch
Locus ID:	54802
UniProt ID:	Q9H3H1



[View online »](#)