

Product datasheet for **TP725237**

GTPBP3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human GTPBP3(N-His-SUMO)
Species:	Human
Expression cDNA Clone or AA Sequence:	H60-L155
Tag:	N-His-SUMO
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human GTPBP3 is produced by E.coli expression system and the target gene encoding H60-L155 is expressed with a His&SUMO tag at the N-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	84705
UniProt ID:	Q969Y2
Synonyms:	tRNA modification GTPase GTPBP3, mitochondrial; GTP-binding protein 3; Mitochondrial GTP-binding protein 1



[View online »](#)