

Product datasheet for **TP725240**

COQ6 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human COQ6 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	N85-E157
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human COQ6 Protein is produced by E.coli expression system and the target gene encoding N85-E157 is expressed.
Stability:	12 months from date of despatch
Locus ID:	51004
UniProt ID:	Q9Y2Z9



[View online »](#)