

Product datasheet for **TP725242**

COQ4 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human COQ4 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	E101-W184
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human COQ4 Protein is produced by E.coli expression system and the target gene encoding E101-W184 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	51117
UniProt ID:	Q9Y3A0
Synonyms:	Ubiquinone biosynthesis protein COQ4 homolog, mitochondrial; Coenzyme Q biosynthesis protein 4 homolog



[View online »](#)