

Product datasheet for **TP725243**

UBE2E2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Ubiquitin-conjugating enzyme E2 E2/UBE2E2
Species:	Human
Expression cDNA Clone or AA Sequence:	S2-C130
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Ubiquitin-conjugating enzyme E2 E2/UBE2E2 is produced by E.coli expression system and the target gene encoding S2-C130 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	7325
UniProt ID:	Q96LR5
Synonyms:	E2 ubiquitin-conjugating enzyme E2; Ubch8; Ubiquitin carrier protein E2; Ubiquitin-protein ligase E2
Protein Pathways:	Ubiquitin mediated proteolysis



[View online »](#)