

## Product datasheet for **TP725246**

### MYO5A Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Unconventional myosin-Va/MYO5A
Species:	Human
Expression cDNA Clone or AA Sequence:	T951-L1150
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Unconventional myosin-Va/MYO5A is produced by E.coli expression system and the target gene encoding T951-L1150 is expressed.
Stability:	12 months from date of despatch
Locus ID:	4644
UniProt ID:	<a href="#">Q9Y4I1</a>



[View online »](#)