

Product datasheet for **TP725255**

C1orf149 (MEAF6) Human Recombinant Protein

Product data:

| | |
|--|--|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant Human Chromatin modification-related protein MEAF6(C-His) |
| Species: | Human |
| Expression cDNA Clone or AA Sequence: | A2-Y191 |
| Tag: | C-His |
| Buffer: | Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4 |
| Note: | MADE TO ORDER Recombinant Human Chromatin modification-related protein MEAF6 is produced by E.coli expression system and the target gene encoding A2-Y191 is expressed with a His tag at the C-terminus. |
| Storage: | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles. |
| Stability: | 12 months from date of despatch |
| Locus ID: | 64769 |
| UniProt ID: | Q9HAF1 |
| Synonyms: | MYST/Esa1-associated factor 6; Esa1-associated factor 6 homolog; Protein EAF6 homolog; hEAF6; Sarcoma antigen NY-SAR-91; EAF6 |



[View online »](#)