

## Product datasheet for **TP725262**

### TAF1L Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human TAF1L Protein(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	L1517-E1649
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human TAF1L Protein is produced by E.coli expression system and the target gene encoding L1517-E1649 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	138474
UniProt ID:	<a href="#">Q8IZX4</a>
Synonyms:	Transcription initiation factor TFIID subunit 1-like; TAF(II)210; TBP-associated factor 1-like; TBP-associated factor 210 kDa; Transcription initiation factor TFIID 210 kDa subunit
Protein Families:	Protein Kinase, Transcription Factors
Protein Pathways:	Basal transcription factors



[View online »](#)