

Product datasheet for **TP725264**

SMC1L2 (SMC1B) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human SMC-1B(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	A855-L982
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human SMC-1B is produced by E.coli expression system and the target gene encoding A855-L982 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	27127
UniProt ID:	Q8NDV3
Synonyms:	Structural maintenance of chromosomes protein 1B;SMC protein 1B; SMC-1-beta; SMC1L2
Protein Pathways:	Cell cycle, Oocyte meiosis



[View online »](#)