

## Product datasheet for **TP725290**

### AARSD1 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Alanyl-tRNA editing protein Aarsd1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	L110-V311
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Alanyl-tRNA editing protein Aarsd1 is produced by E.coli expression system and the target gene encoding L110-V311 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	80755
UniProt ID:	<a href="#">Q9BTE6</a>



[View online »](#)