

Product datasheet for **TP725293**

Gkn1 Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Mouse Gastrokine-1(C-His)
Species:	Mouse
Expression cDNA Clone or AA Sequence:	V26-N137
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Mouse Gastrokine-1 is produced by E.coli expression system and the target gene encoding V26-N137 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	66283
UniProt ID:	Q9CR36



[View online »](#)