

Product datasheet for **TP725294**

EPS8L2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human EPS8-like protein 2(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	S366-M493
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human EPS8-like protein 2 is produced by E.coli expression system and the target gene encoding S366-M493 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	64787
UniProt ID:	Q9H6S3
Synonyms:	Epidermal growth factor receptor kinase substrate 8-like protein 2; Epidermal growth factor receptor pathway substrate 8-related protein 2; EPS8-related protein 2;EPS8L2; EPS8R2
Protein Families:	Druggable Genome



[View online »](#)