

Product datasheet for **TP725298**

SIRT5 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human SIRT5 Protein(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	P95-C242
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human SIRT5 Protein is produced by E.coli expression system and the target gene encoding P95-C242 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	23408
UniProt ID:	Q9NXA8
Synonyms:	NAD-dependent protein deacylase sirtuin-5, mitochondrial; Regulatory protein SIR2 homolog 5; SIR2-like protein 5; SIR2L5
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)