

## Product datasheet for **TP725300**

### **RING1 Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Human E3 ubiquitin-protein ligase RING1(C-His)
<b>Species:</b>	Human
<b>Expression cDNA Clone or AA Sequence:</b>	R78-I363
<b>Tag:</b>	C-His
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
<b>Note:</b>	MADE TO ORDER   Recombinant Human E3 ubiquitin-protein ligase RING1 is produced by E.coli expression system and the target gene encoding R78-I363 is expressed with a His tag at the C-terminus.
<b>Storage:</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
<b>Stability:</b>	12 months from date of despatch
<b>Locus ID:</b>	6015
<b>UniProt ID:</b>	<a href="#">Q06587</a>
<b>Synonyms:</b>	Polycomb complex protein RING1; RING finger protein 1; RING-type E3 ubiquitin transferase RING1; Really interesting new gene 1 protein; RNF1
<b>Protein Families:</b>	Druggable Genome, Transcription Factors



[View online »](#)