

Product datasheet for **TP725304**

PCLKC (CDHR2) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Cadherin-related family member 2/CDHR2(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	D763-P903
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Cadherin-related family member 2/CDHR2 is produced by E.coli expression system and the target gene encoding D763-P903 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	54825
UniProt ID:	Q9BYE9
Synonyms:	Protocadherin LKC; Protocadherin-24
Protein Families:	Transmembrane



[View online »](#)