

Product datasheet for **TP725305**

NT5C Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human NT5C Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-E201
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human NT5C Protein is produced by E.coli expression system and the target gene encoding M1-E201 is expressed.
Stability:	12 months from date of despatch
Locus ID:	30833
UniProt ID:	Q8TCD5



[View online »](#)