

## Product datasheet for **TP725307**

### SLC9A3R2 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human NHERF-2
Species:	Human
Expression cDNA Clone or AA Sequence:	L56-F337
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human NHERF-2 is produced by E.coli expression system and the target gene encoding L56-F337 is expressed.
Stability:	12 months from date of despatch
Locus ID:	9351
UniProt ID:	<a href="#">Q15599</a>



[View online »](#)