

Product datasheet for **TP725309**

Cytoglobin (CYGB) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Cytoglobin/CYGB(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-P190
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Cytoglobin/CYGB is produced by E.coli expression system and the target gene encoding M1-P190 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	114757
UniProt ID:	<u>Q8WWM9</u>



[View online »](#)