

## Product datasheet for **TP725310**

### Protein phosphatase inhibitor 1 (PPP1R1A) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Protein phosphatase 1 regulatory subunit 1A(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-V171
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Protein phosphatase 1 regulatory subunit 1A is produced by E.coli expression system and the target gene encoding M1-V171 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	5502
UniProt ID:	<a href="#">Q13522</a>



[View online »](#)