

Product datasheet for **TP725314**

MOBK1B (MOB1A) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human MOB kinase activator 1A/MOB1A(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	E33-R216
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human MOB kinase activator 1A/MOB1A is produced by E.coli expression system and the target gene encoding E33-R216 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	55233
UniProt ID:	Q9H8S9



[View online »](#)