

## Product datasheet for **TP725325**

### LSM4 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human U6 snRNA-associated Sm-like protein LSM4(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-Q139
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human U6 snRNA-associated Sm-like protein LSM4 is produced by E.coli expression system and the target gene encoding M1-Q139 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	25804
UniProt ID:	<a href="#">Q9Y4Z0</a>
Synonyms:	Glycine-rich protein; GRP
Protein Families:	Stem cell - Pluripotency
Protein Pathways:	RNA degradation, Spliceosome



[View online »](#)