

Product datasheet for **TP725364**

UBF1 (UBTF) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Nucleolar transcription factor 1/UBTF(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	K200-L440
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Nucleolar transcription factor 1/UBTF is produced by E.coli expression system and the target gene encoding K200-L440 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	7343
UniProt ID:	P17480
Synonyms:	Autoantigen NOR-90; Upstream-binding factor 1; UBF-1
Protein Families:	Transcription Factors



[View online »](#)