

Product datasheet for **TP725366**

H10 (H1F0) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Histone H1.0(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	T2-K194
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Histone H1.0 is produced by E.coli expression system and the target gene encoding T2-K194 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	3005
UniProt ID:	P07305
Synonyms:	Histone H1'; Histone H1(0); Histone H1.0, N-terminally processed; H1-0; H1F0; H1FV
Protein Families:	Druggable Genome



[View online »](#)