

## Product datasheet for **TP725378**

### SALL1 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Sal-like protein 1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	P558-F670
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Sal-like protein 1 is produced by E.coli expression system and the target gene encoding P558-F670 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	6299
UniProt ID:	<a href="#">Q9NSC2</a>
Synonyms:	Spalt-like transcription factor 1; Zinc finger protein 794; Zinc finger protein SALL1; Zinc finger protein Spalt-1; HSal1; Sal-1; SALL1; ZNF794
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)