

Product datasheet for **TP725381**

OVGP1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Oviduct-specific glycoprotein/OVGP1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	S481-A678
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Oviduct-specific glycoprotein/OVGP1 is produced by E.coli expression system and the target gene encoding S481-A678 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	5016
UniProt ID:	Q12889
Synonyms:	Estrogen-dependent oviduct protein; Mucin-9; Oviductal glycoprotein; Oviductin;MUC9; OGP
Protein Families:	Secreted Protein



[View online »](#)