

Product datasheet for **TP725383**

KRT25 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Keratin, type I cytoskeletal 25/ KRT25(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	E79-C394
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Keratin, type I cytoskeletal 25/ KRT25 is produced by E.coli expression system and the target gene encoding E79-C394 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	147183
UniProt ID:	Q7Z3Z0
Synonyms:	Cytokeratin-25; CK-25; Keratin-25; K25; Keratin-25A; K25A; Type I inner root sheath-specific keratin-K25irs1; KRT25A



[View online »](#)