

## Product datasheet for **TP725394**

### SFRS8 (SFSWAP) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human SFSWAP Protein(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	S170-G280
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human SFSWAP Protein is produced by E.coli expression system and the target gene encoding S170-G280 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	6433
UniProt ID:	<a href="#">Q12872</a>
Synonyms:	Splicing factor, suppressor of white-apricot homolog; Splicing factor, arginine/serine-rich 8; Suppressor of white apricot protein homolog; SFRS8; SWAP
Protein Families:	Transcription Factors



[View online »](#)