

Product datasheet for **TP725403**

HOX11 (TLX1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human T-cell leukemia homeobox protein 1/TLX1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	K201-E330
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human T-cell leukemia homeobox protein 1/TLX1 is produced by E.coli expression system and the target gene encoding K201-E330 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	3195
UniProt ID:	P31314
Synonyms:	Homeobox protein Hox-11; Proto-oncogene TCL-3; T-cell leukemia/lymphoma protein 3



[View online »](#)