

## Product datasheet for **TP725410**

### NAPE PLD (NAPEPLD) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human NAPE-PLD(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	D2-L125
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human NAPE-PLD is produced by E.coli expression system and the target gene encoding D2-L125 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	222236
UniProt ID:	<a href="#">Q6IQ20</a>
Synonyms:	N-acyl-phosphatidylethanolamine-hydrolyzing phospholipase D; N-acyl phosphatidylethanolamine phospholipase D; NAPE-hydrolyzing phospholipase D



[View online »](#)