

Product datasheet for **TP725418**

CtIP (RBBP8) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human DNA endonuclease RBBP8(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	K604-T897
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human DNA endonuclease RBBP8 is produced by E.coli expression system and the target gene encoding K604-T897 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	5932
UniProt ID:	Q99708
Synonyms:	CtBP-interacting protein; CtIP; Retinoblastoma-binding protein 8; RBBP-8; Retinoblastoma-interacting protein and myosin-like; RIM; Sporulation in the absence of SPO11 protein 2 homolog; SAE2



[View online »](#)