

## Product datasheet for **TP725427**

### SMARCE1 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human SMARCE1 Protein(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	D90-E411
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human SMARCE1 Protein is produced by E.coli expression system and the target gene encoding D90-E411 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	6605
UniProt ID:	<a href="#">Q969G3</a>
Synonyms:	SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily E member 1; BRG1-associated factor 57; BAF57
Protein Families:	Transcription Factors



[View online »](#)