

Product datasheet for **TP725431**

WDHD1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human WDHD1 Protein(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	E830-E1129
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human WDHD1 Protein is produced by E.coli expression system and the target gene encoding E830-E1129 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	11169
UniProt ID:	O75717
Synonyms:	WD repeat and HMG-box DNA-binding protein 1; Acidic nucleoplasmic DNA-binding protein 1; And-1
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)