

Product datasheet for **TP725432**

BRD7 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Bromodomain-containing protein 7/BRD7(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	G2-T350
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Bromodomain-containing protein 7/BRD7 is produced by E.coli expression system and the target gene encoding G2-T350 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	29117
UniProt ID:	Q9NPI1



[View online »](#)