

Product datasheet for **TP725437**

KHDRBS3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human KHDRBS3 Protein(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	A11-D300
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human KHDRBS3 Protein is produced by E.coli expression system and the target gene encoding A11-D300 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10656
UniProt ID:	O75525
Synonyms:	KH domain-containing, RNA-binding, signal transduction-associated protein 3; RNA-binding protein T-Star; Sam68-like mammalian protein 2; Sam68-like phosphotyrosine protein; SALP; SLM2



[View online »](#)