

Product datasheet for **TP725440**

PAK1 interacting protein 1 (PAK1IP1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human PAK1IP1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	N51-T390
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human PAK1IP1 is produced by E.coli expression system and the target gene encoding N51-T390 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	55003
UniProt ID:	Q9NWT1
Synonyms:	p21-activated protein kinase-interacting protein 1; PAK/PLC-interacting protein 1; hPIP1; PAK1-interacting protein 1; WD repeat-containing protein 84



[View online »](#)