

Product datasheet for **TP725441**

LACRT Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human LACRT/Lacritin (N-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	Ala19-Ala138
Tag:	N-6His
Buffer:	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4
Note:	Recombinant Human Lacritin is produced by our E.coli expression system and the target gene encoding Ala19-Ala138 is expressed with a 6His tag at the N-terminus.
Stability:	12 months from date of despatch
Locus ID:	90070
UniProt ID:	Q9GZZ8
Summary:	Human LACRT is primarily localized to secretory granules and fluid. This protein is highly expressed in the lacrimal glands. This protein promotes ductal cell proliferation and stimulates signaling through tyrosine phosphorylation.



[View online »](#)