

Product datasheet for **TP725442**

MESP1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Mesoderm posterior protein 1/MESP1(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	P21-L250
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Mesoderm posterior protein 1/MESP1 is produced by E.coli expression system and the target gene encoding P21-L250 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	55897
UniProt ID:	Q9BRJ9
Synonyms:	Class C basic helix-loop-helix protein 5; bHLHc5
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)