

## Product datasheet for **TP725443**

### **SPT3 (SUPT3H) Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Human Transcription initiation protein SPT3 homolog(C-His)
<b>Species:</b>	Human
<b>Expression cDNA Clone or AA Sequence:</b>	N2-A230
<b>Tag:</b>	C-His
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
<b>Note:</b>	MADE TO ORDER   Recombinant Human Transcription initiation protein SPT3 homolog is produced by E.coli expression system and the target gene encoding N2-A230 is expressed with a His tag at the C-terminus.
<b>Stability:</b>	12 months from date of despatch
<b>Locus ID:</b>	8464
<b>UniProt ID:</b>	<a href="#">O75486</a>



[View online »](#)