

Product datasheet for **TP725444**

ARNT2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human ARNT protein 2(C-His)
Species:	Human
Expression cDNA Clone or AA Sequence:	Q421-Q650
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human ARNT protein 2 is produced by E.coli expression system and the target gene encoding Q421-Q650 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	9915
UniProt ID:	Q9HBZ2
Synonyms:	Aryl hydrocarbon receptor nuclear translocator 2; Class E basic helix-loop-helix protein 1; BHLHE1
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Pathways in cancer, Renal cell carcinoma



[View online »](#)