

Product datasheet for **TP725448**

DHX38 Human Recombinant Protein

Product data:

| | |
|--|--|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant Human DHX38 Protein(C-His) |
| Species: | Human |
| Expression cDNA Clone or AA Sequence: | D19-D280 |
| Tag: | C-His |
| Buffer: | Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4 |
| Note: | MADE TO ORDER Recombinant Human DHX38 Protein is produced by E.coli expression system and the target gene encoding D19-D280 is expressed with a His tag at the C-terminus. |
| Stability: | 12 months from date of despatch |
| Locus ID: | 9785 |
| UniProt ID: | <u>Q92620</u> |



[View online »](#)